

LYMPHA PRESS-IVE!

LYMPHA PRESS-IVE INDEED! Yes, it's quite impressive what many of the Medical Industry's experts are saying about the Lympha Press Pump's multi-compartmental, calibrated, gradient pressure system specifically designed for the effective treatment of lymphatic and venous disorders.



Primary Lymphedema (Model 103AT)



Secondary Lymphedema (Model 201A)

A Lympha Press 201A is recommended for those patients whose ulcers have not healed as well as the small group of patients with the following handicap(s):

- Hip replacement
- arthritis of the hands
- any patients that cannot dress themselves
- wheelchair confinement
- obesity

These patients have lost the ability to apply the compression gradient stockings that doctors have prescribed for them. Consequently, these patients will always need a Lympha Press Pump for which we recommend our 201A Model. We will ensure that those in need will receive the Lympha Press Pump through our network of certified technicians. Our Lympha Press Pump instructional product videos are available for your viewing anytime – you can order one by contacting us as follows:

LYMPHA PRESS -PRESS!

The protocol with the use of Lympha Press is available for the following diagnoses:

- Pediatric Lymphedema
- Lymphedema of the genitals
- Lymphedema (primary/secondary)
- Excisional Surgery
- Venous Ulcer

The age range of our protocol is anywhere from newborns at 5 weeks of age to seniors at 96 years of age.

The Lympha Press Pump offers patients many features and benefits. Although they are too numerous to mention them all here – the following are a few “key highlights” of this product:

1 Movement of fluid in a distal to proximal direction.

The unique overlapping of the air compartments within the sleeve guarantees a calibrated gradient pressure and prevents any gaps between the functioning compartments creating a “milking” effect on the affected limb which produces nearly the same quality of lymphatic fluid as is present within a healthy limb.

2 Multi-compartmental pressure system.

This system supplies pressure from the distal area of the sleeve and progresses proximally through the overlapping cells thereby providing a calibrated gradient “milking” pattern.

3 Short 30-second cycles.

These brief cycles allow for comfortable, high-pressure treatment to encourage maximum patient compliance and comfort.

4 Optional variable pressure controls.

Provides control over calibrated gradient pressure from 30-200mm Hg.

Regardless of diagnosis – early intervention and long-term maintenance using the Lympha Press is recognized by many of today’s leading physicians as an effective treatment to be used towards bringing the patient’s limb back as close as possible to its natural condition. The Lympha Press is an important aid to the physician before and after the surgical removal of scar tissue. It can reduce the hospitalization period and post-operative infection rate while minimizing pain.